

# CYTOSIAL

DEEP Wrinkles MEDIUM VOLUME

## SKIN REJUVENATION

Improvement of skin hydration, firmness, tone and elasticity	Face			
	Decolletage		✓	
	Hands		✓	

## EFFECTIVE TREATMENT OF LINES AND FOLDS

Mild and superficial lines	Crow's feet		✓	
	Oral commissure	✓	✓	
	Forehead lines		✓	
	Glabellar lines		✓	
Moderate folds	Nasolabial folds	✓	✓	
	Marionette lines	✓	✓	
	Mental crease		✓	
Deep and severe folds	Marionette lines		✓	
	Mental crease	✓	✓	

## LIP AUGMENTATION

	Lip contouring, vermillion border	✓	✓	
	Volumisation of upper and lower lip	✓		✓
	Accentuation of Cupid's bow		✓	

## VOLUME ENHANCEMENT

	Temporal hollows		✓	✓
	Medial cheek/ Lateral cheek		✓	✓
	Submalar cheek hollows		✓	
	Mandibular angle		✓	✓
	Earlobes		✓	✓

Recommended Treatment Pattern

BIOXIS











info@amittismed.com  
amittismed.com  
amittismedco

(021) 22 90 90 54  
0902 285 38 80

# Cytosial (Product Range)

**CYTOSIAL**   
The filler you love

Monophasic polydensified products

Product	HA concentration/ cross - linking level	Needle size	Volume	Indications	level of injection
<p><b>VOLUME</b></p> 	HA : 25mg/ ml	27G <sup>1/2</sup>	1.1ml	<p>Restore facial volume:</p> <ul style="list-style-type: none"> <li>•cheeks and cheekbones</li> <li>•chin</li> <li>•facial volume loss</li> <li>•correction of the shape of the nose</li> <li>•line and angle of the jaw</li> <li>•contour and volume of the beard</li> <li>•cheekbones</li> </ul> 	<p>-Deep dermis or subcutaneous</p> 
<p><b>DEEP Wrinkles</b></p> 	HA : 25mg/ ml	27G <sup>1/2</sup>	1.1ml	<p>Correct deep lines:</p> <ul style="list-style-type: none"> <li>•deep nasolabial fold</li> <li>•deep marionette fold</li> <li>•deep facial wrinkles</li> <li>•the middle part of the face, valley of tears, cheeks</li> <li>•around the temple and eyebrows</li> <li>•nose and lip</li> <li>•beard</li> </ul> 	<p>-Deep dermis</p> 
<p><b>MEDIUM</b></p> 	HA : 25mg/ ml	30G <sup>1/2</sup>	1.1ml	<p>Is indicated for:</p> <ul style="list-style-type: none"> <li>•the filling of moderate cutaneous depressions as perioral or commissural areas</li> <li>•mouth area (contour, volume)</li> <li>•facial volume loss</li> <li>•modeling of the ear lobes</li> <li>•lip contour and perioral line</li> </ul> 	<p>-Superficial to Mid dermis</p> 